

## Discovery of Variants & Polymorphisms

### Applications

Genetic variation forms the foundation of phenotypic diversity, trait inheritance, and adaptive evolution. Theracues' Variant and Polymorphism Discovery service enables comprehensive identification of single nucleotide polymorphisms (SNPs), simple sequence repeats (SSRs), and repeat elements using high-throughput sequencing technologies, supporting applications in genomics, breeding, population genetics, and molecular marker development.

#### Pros

- Genome-wide detection of SNPs, SSRs, and repeats
- High sensitivity and accuracy
- Supports population and comparative studies
- Enables marker-assisted selection
- Applicable across multiple species

#### Notes

- Requires adequate sequencing depth
- Reference genome improves accuracy
- Complex genomes may need higher coverage
- Validation recommended for functional use

### Input



Genomic DNA

DNA Quantity: ≥200 ng  
(recommended)

### Sample Types



Blood/  
PBMC



Tissue



Plant tissues  
(leaf, root,  
stem)




Bacterial  
Cultures

**Platform** Illumina/ ONT Long-read sequencing platforms

### Process Map



Sample  
DNA  
Extraction



DNA  
QC



Library  
Prep



NGS Data  
Acquisition



Data QC



Read  
Alignment



Variant calling  
(SNPs, InDels, SSRs,  
repeats)



Variant  
annotation

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.

## Discovery of Variants & Polymorphisms

### Applications

Rapid Molecular Diversity Analysis enables quick and reliable assessment of genetic variation within and between populations. By leveraging molecular markers and high-throughput sequencing or targeted genotyping approaches, this service supports diversity studies, germplasm characterization, breeding programs, and evolutionary research. This service focuses on evaluating genetic diversity, population structure, and relatedness using molecular marker data. It enables rapid comparison of multiple samples to identify genetic similarities and differences, supporting applications in breeding selection, conservation genetics, and population-level studies. Analyses can be performed using SNPs, SSRs, or other polymorphic markers, depending on study goals and turnaround requirements.

#### Pros

- Rapid turnaround time
- Efficient population-level analysis
- High reproducibility
- Supports breeding and conservation studies
- Scalable across sample sizes

#### Notes

- Marker choice impacts resolution
- Reference data improves interpretation
- Limited functional inference
- Best suited for comparative analysis

### Input



Genomic DNA

DNA Quantity:  $\geq 100$  ng  
(recommended)

### Sample Types



Blood/  
PBMC



Tissue



Plant tissues  
(leaf, root,  
stem)




Bacterial  
Cultures

**Platform** Illumina sequencing platforms

### Process Map



Sample  
DNA  
Extraction



DNA  
QC



Library  
Prep



NGS Data  
Acquisition



Data QC



Read  
Alignment



Variant calling  
(SNPs, InDels, SSRs,  
repeats)



Variant  
annotation



Diversity and  
population  
analysis

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.

## Cell line characterization - clonality analysis

### Applications

Cell line authentication is essential to ensure research validity, reproducibility, and data integrity. 'Cell Line Typing - Amplicon' service provides accurate genetic profiling of cell lines using targeted amplicon sequencing, enabling reliable identification, authentication, and detection of cross-contamination. This service uses amplicon-based sequencing to analyze specific genetic loci for cell line identification and verification. It supports routine authentication, quality control, and validation of cell lines used in research and translational studies. The workflow enables comparison against reference profiles to confirm identity and detect misidentification or contamination, ensuring experimental reliability.

#### Pros

- High accuracy and reproducibility
- Rapid and cost-effective
- Detects cross-contamination
- Suitable for routine QC
- Low DNA input required

#### Notes

- Limited to targeted loci
- Reference profiles required
- Not genome-wide
- Periodic testing recommended

### Input



Genomic DNA

DNA Quantity:  $\geq 100$  ng  
(recommended)

### Sample Types



Human  
Cell-lines



Animal  
Cell-lines



Fresh or  
frozen cell  
pellets

**Platform** Illumina Sequencing Platforms

### Process Map

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.

## Hybridoma sequencing

### Applications

This service focuses on sequencing antibody variable domains (VH and VL) from hybridoma-derived samples. By combining targeted amplification with next-generation sequencing, it enables reliable recovery of functional antibody sequences for cloning, expression, humanization, and intellectual property protection. Hybridoma sequencing ensures long-term preservation of antibody sequence information independent of cell line stability.

#### Pros

- Accurate VH/VL sequence recovery
- Supports antibody cloning
- Preserves antibody IP
- High sensitivity and specificity
- Suitable for R&D workflows

#### Notes

- Requires high-quality RNA/DNA
- Targets variable regions only
- Mixed clones may complicate analysis
- Confirmation recommended

### Input



Total RNA or Genomic DNA

- RNA:  $\geq 500$  ng
- DNA:  $\geq 200$  ng

### Sample Types



Hybridoma  
cell pellets



Cultured  
hybridoma  
cells

**Platform** Illumina/ ONT Long-read sequencing platforms

### Process Map

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.

## Bulk Immune Repertoire Sequencing (B-cell / T-cell)

### Applications

This service profiles B-cell receptors (BCRs) and T-cell receptors (TCRs) using long-read sequencing technologies to recover full-length receptor sequences in bulk samples. It enables detailed analysis of:

- Immune diversity and clonotype composition
- Clonal expansion and dominance
- Isotype and chain pairing (BCR)
- Immune responses in disease, vaccination, and therapy

Long-read sequencing eliminates assembly ambiguity associated with short reads, ensuring accurate repertoire reconstruction.

#### Pros

- Full-length BCR/TCR sequences
- Accurate clonotype resolution
- Reduced assembly bias
- Captures V(D)J and CDR regions
- Suitable for complex repertoires

#### Notes


- Requires high-quality RNA
- Bulk-level resolution only
- Higher input recommended
- Data interpretation is specialized

### Input



Total RNA Input Amount:  $\geq 1 \mu\text{g}$

RNA Integrity (RIN)



### Sample Types



PBMC's



Whole  
Blood



Sorted B-cell  
populations



Lymphoid  
tissues

**Platform** ONT Long-read sequencing platforms

### Process Map

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.

## BCR Sequencing

### Applications

This service focuses on sequencing immunoglobulin heavy and light chain regions to analyze B-cell repertoires in bulk samples. It enables assessment of:

- B-cell diversity and clonotype distribution
- Clonal expansion and dominance
- V(D)J gene usage and CDR characteristics
- Immune responses in infection, vaccination, and disease

The workflow supports both short-read and long-read sequencing strategies, depending on resolution requirements.

#### Pros

- High-resolution BCR profiling
- Sensitive detection of clonal expansion
- Supports immunology and vaccine studies
- Flexible sequencing strategies
- Reproducible and scalable

#### Notes


- Bulk-level analysis only
- RNA quality is critical
- Chain pairing limited with short reads
- Specialized interpretation required

### Input



Total RNA Input Amount:  $\geq 500$  ng

RNA Integrity (RIN)



### Sample Types



PBMC's



Whole  
Blood



Sorted B-cell  
populations



Lymphoid  
tissues

**Platform** ONT Long-read sequencing platforms

### Process Map

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.

## TCR Sequencing

### Applications

This service sequences TCR  $\alpha$ ,  $\beta$  (and optionally  $\gamma$ ,  $\delta$ ) chains to analyze T-cell repertoires in bulk samples. It enables assessment of:

- T-cell diversity and clonotype composition
- Clonal expansion and immune dominance
- V(D)J gene usage and CDR3 diversity
- Immune monitoring in disease, vaccination, and therapy

The workflow supports short-read or long-read sequencing, depending on resolution and pairing requirements.

#### Pros

- High-resolution TCR profiling
- Sensitive clonotype detection
- Supports immuno-oncology studies
- Scalable across cohorts
- Reproducible workflows

#### Notes


- Bulk-level analysis only
- RNA quality is critical
- Chain pairing limited with short reads
- Expert interpretation required

### Input



Total RNA Input Amount:  $\geq 1 \mu\text{g}$

RNA Integrity (RIN)



### Sample Types



PBMC's



Whole  
Blood



Sorted B-cell  
populations



Lymphoid  
tissues

**Platform** ONT Long-read sequencing platforms

### Process Map

### Standard Deliverables

[Download Sample Report](#)

### References

- Ellegren, H. (2004). Microsatellites: simple sequences with complex evolution. *Nature Reviews Genetics*.
- Nielsen, R. et al. (2011). Genotype and SNP calling from next-generation sequencing data. *Nature Reviews Genetics*.